











Soft-Release-O™

Infusion Set



If self-injecting only.

click!

Soft-Release-O™

en	Instructions for use	10
ar	إرشادات الاستخدام	14
cs	Návod k použití	18
da	Brugsanvisning	22
de	Gebrauchsanweisung	26
el	Οδηγίες χρήσης	30
es	Instrucciones de uso	34
et	Kasutusjuhend	38
fi	Käyttöohje	42
fr	Mode d'emploi	46
he	הוראות שימוש	50
hr	Upute za upotrebu	54
hu	Használati útmutató	58
it	Istruzioni d'uso	62
ja	使用説明	66
ko	사용 설명서	70
lt	Naudojimo instrukcija	74
lv	Lietošanas instrukcija	78
nl	Gebruiksaanwijzing	82
no	Bruksinstruksjoner	86
pl	Instrukcja obsługi	90
pt	Instruções de utilização	94
ro	Instrucțiuni de utilizare	98
ru	Инструкция по использованию	102
sk	Návod na použitie	106
sl	Navodila za uporabo	110
sv	Bruksanvisning	114
tr	Kullanım talimatları	118
zh	使用说明	122

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English

Indications for Use

Soft-Release-O™ Infusion Sets are intended for the subcutaneous delivery of medication, such as insulin, from an external infusion pump.

Contraindications

Soft-Release-O™ Infusion Sets are not intended for intravenous (IV) infusion or the infusion of blood or blood products.

Warnings

- ⚠ **Follow the instructions for use** of Soft-Release-O™ Infusion Set and your infusion pump. Use the Soft-Release-O™ Infusion Set with your infusion pump only after consultation with your physician and/or health care professional.
- ⚠ **Avoid infections.** Do not use the Soft-Release-O™ Infusion Set if the individual packaging is damaged, has already been opened, or if the cannula cap is missing. Do not use the Soft-Release-O™ Infusion Sets after the expiration date on the packaging. If these instructions are not followed, there is a risk of infection. Seek immediate medical attention if you have any health concern.
- ⚠ **Single use only!** Do not re-use, re-insert and/or re-sterilize the Soft-Release-O™ Infusion Set. Re-use and/or re-insertion may lead to infection, site irritation, inaccurate medication delivery and/or may cause damage to the cannula. Re-sterilization may affect the integrity of the Soft-Release-O™ Infusion Set. Seek immediate medical attention if you have any health concern.
- ⚠ **Do not use Soft-Release-O™ Infusion Set for longer than 72 hours.** Longer use may lead to infection, site irritation and/or inaccurate medication delivery. Follow your physician's instructions for the amount of time the Soft-Release-O™ Infusion Set can be used. Change your infusion site if it becomes itchy, red, inflamed or if you have an unexplained high blood sugar reading.
- ⚠ **Ensure correct delivery of medication.** Before insertion, check that the Soft-Release-O™ Infusion Set allows free flow of medication, following the priming instructions for your infusion pump. This will prevent incorrect dosage caused by air in the fluid path or a clogged needle. Use a new Soft-Release-O™ Infusion Set if priming is not successful.
Never prime the Soft-Release-O™ Infusion Set or attempt to free a clogged line while the Soft-Release-O™ Infusion Set is inserted. You may accidentally inject too much medication. Do not re-insert the introducer needle into the Soft-Release-O™ Infusion Set. Re-insertion may cause tearing or puncturing of the soft cannula, which may result in unpredictable medication flow. Replace the Soft-Release-O™ Infusion Set if the tape becomes loose, if the soft cannula becomes kinked during insertion or if

Soft-Release-O™

the cannula becomes fully or partially dislodged from the skin. Check frequently to make sure the soft cannula remains firmly in place. It must always be completely inserted to receive the full amount of medication.

Do not change the Soft-Release-O™ Infusion Set just before bedtime. A blood glucose test has to be performed 1-2 hours after insertion to ensure adequate delivery of medication.

Develop a conventional insulin injection or medication replacement plan with your physician if any part of the infusion system fails or system supplies are depleted.

If these instructions are not followed there is a risk of incorrect delivery of medication. Seek immediate medical attention if you have any health concern.

- ⚠ **Avoid bending or breaking** the introducer needle. A broken needle can remain stuck in the body or remain completely under the skin. If this occurs seek immediate medical attention.
- ⚠ For Healthcare Professionals: **Avoid handling, recapping or destroying** the introducer needle. Handling, recapping or destroying the introducer needle can lead to accidental pricking which can transmit serious infections. Professional health care personnel should follow the specific directions for health care personnel including, without limitation, OSHA (Occupational Safety and Health Administration), CDC (Center for Disease Control and Prevention), or local policy when using infusion set, including directions for recapping and disposing. To obtain the specific directions, consult a physician. Seek immediate medical attention if you did prick yourself.

Caution (please note)

- Soft-Release-O™ Infusion Sets are designed to work with the DANA DIABECARE insulin pump only. Infusion pump and/or pump reservoir manufacturer's guide must be followed.
- Change the infusion site regularly in accordance with your physician's and/or health care professional's advice. If this is not done you risk insulin resistance, inflammation of the infusion site or lipohypertrophy.
- Consult your physician who will recommend the correct cannula length. It is important to match the cannula length to your body physique and the selected infusion site in order to ensure accurate medication delivery.
- Do not put disinfectants, perfumes, deodorants or insect repellent on the Soft-Release-O™ Infusion Set. Protect the product from excessive humidity, sunlight and heat. This may affect the integrity of the Soft-Release-O™ Infusion Set. Store Soft-Release-O™ Infusion Sets at room temperature.

English

Disposal

You must follow your physician's and/or health care professional's advice and local government regulations for disposal of all used needles/sharps.

Figure 1-2

Wash hands and practice aseptic technique. Prepare the infusion site by removing excess hair. Using a 70% isopropyl alcohol pad, wipe the site in a circular outward motion. Make sure site is completely dry before proceeding. Failure to use 70% isopropyl alcohol pad or use of other than 70% isopropyl alcohol pads may result in premature tape failure and interruption of medication delivery.

Figure 3

After connection to the infusion pump make sure all air is removed by priming the tubing assembly according to your pump manufacturer's priming instructions.

Figure 4-5

Remove tape backing and be careful not to touch adhesive. Carefully remove cannula cap.

Figure 6-7

Stabilize the infusion site and insert Soft-Release-O™ cannula with an angle of 90°. For maximum adhesion, press and hold tape to skin for 15 seconds to allow tape to reach body temperature.

Figure 8-10

Hold tape against skin with two fingers of one hand. Gently squeeze tabs and remove introducer needle carefully with the other hand. Dispose the introducer needle into a safe container unit immediately after injection to avoid accidentally pricking. Self injecting users may use the blue needle cover provided to temporarily cover the introducer needle after removal.

Figure 11-13

Swab septum with a 70% isopropyl alcohol pad and the small plastic cannula on the bottom of the tube cap with the corner of a 70% isopropyl alcohol pad. Proper swabbing of septum & tube cap cannula reduces germ ingress into the fluid path. Attach tube cap to the infusion base. One or two clicks should be heard. Rotate the tubing left and right, at least one full turn in each direction, while pulling upward on cap to ensure the cap is fully engaged.

Figure A: Disconnection

To disconnect the tubing, gently squeeze tabs on both sides of the cap and remove it. There is no need to cover the infusion base as the septum provides closure.

Before you reconnect the cap remove all air from the tubing. To reconnect, follow instructions in figures 11-13.

Monitor your blood glucose levels when disconnected and after reconnecting.

Priming volume for the tubing assembly (volumes are approximate):

30 cm (12 inch): 7 units* (0.07 ml)

60 cm (24 inch): 14 units* (0.14 ml)

80 cm (31 inch): 19 units* (0.19 ml)

100 cm (39 inch): 22 units* (0.22 ml)

Priming volume for the cannula assembly (volumes are approximate):

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Ypsomed AG
Sooil Development Co., Ltd:
Sooil, LLC: 5665 Oberlin Drive, Suite 112,