



INNOVATIVE SOLUTIONS
FOR DIABETES CARE

Philosophy.....	04
Company Profile.....	06
Blood Glucose Monitoring System	08
Professional Care – HbA1c.....	11
Insulin Pump	12
Insulin Therapy – Infusion Sets	13
Insulin Therapy – Pen Needles	14
Lancing Device.....	16
Blood Lancets.....	17
Safety Product	18
Skin Care Product	20
Disinfection & Hygienic Disposal	22
Software Solution	23
What we are looking for.....	24
What we offer	25
References.....	26






IME-DC's objective is to achieve best care for people with diabetes by providing innovative solutions.

We produce and market high-quality products to be able to focus on the most important thing in the world:

The Health of People!

COMPANY PROFILE

- Manufacturer of Blood Glucose Meters as well as Blood Glucose Test Strips
- Patented Technologies, Trademark Owner
- Research & Development Team, Design Department, Laboratory, Production, Sales and Service Team
- Specialised team of 200 employees worldwide
- Wide range of established and well-known brands of diabetes care products
- Proprietary brand represented in 35 countries
- High technology and highest quality manufacturer
- Significant market shares in several countries





HEADQUARTERS
GERMANY Bavaria
Hof

BLOOD GLUCOSE MONITORING SYSTEM

IME-DC

IME-DC Blood Glucose Meter

- Measuring Units available in mg/dL or mmol/L
- Large Icons for a better readability
- Easy Handling for a faster measurement
- Meets the System Accuracy according to EN ISO 15197:2015

IME-DC Blood Glucose Test Strips

- GOD Enzyme enables the blood glucose measurement
- Air-Channel-Sip-In Technology facilitates the blood supply
- Available packaging units: 10 pcs. | 25 pcs. | 2x 25 pcs. | 50 pcs. | 2x 50 pcs.

▼ Product video:



Meets the system accuracy according to EN ISO 15197:2015

BLOOD GLUCOSE MONITORING SYSTEM

DEXTRA Blood Glucose Meter

- Measuring Units available in mg/dL or mmol/L
- Colours available in white or dark blue
- Large Icons for a better readability
- Easy Handling for a faster measurement
- **7 PLUS** quality
- Meets the System Accuracy according to EN ISO 15197:2015



DEXTRA Blood Glucose Test Strips

- GDH-FAD Enzyme ensures a more accurate measurement
- Air-Channel-Sip-In Technology facilitates the blood supply
- Available packaging units: 10 pcs. | 25 pcs. | 2x 25 pcs. | 50 pcs. | 2x 50 pcs.

▼ Product video:



Meets the system accuracy according to EN EN ISO 15197:2015

BLOOD GLUCOSE MONITORING SYSTEM

iDia

NEW

iDia Blood Glucose Meter

- Measuring Units available in mg/dL or mmol/L
- Colours available in white or dark grey (in 2020)
- Large Icons for a better readability
- Easy Handling for a faster measurement
- Display Light for reading in the dark
- Latest technology
- Meets the System Accuracy according to EN ISO 15197:2015

iDia Blood Glucose Test Strips

- GDH-FAD Enzyme ensures a more accurate measurement
- Air-Channel-Sip-In Technology facilitates the blood supply
- Available packaging units: 10 pcs. | 25 pcs. | 2x 25 pcs. | 50 pcs. | 2x 50 pcs.



Meets the
system accuracy
according to
EN ISO
15197:2015

Quo-Lab POCT Analyser¹⁾

- Analyser for HbA1c measurement
- Semi-automated analyser
- Easy to use (step by step instructions on display)
- Results in just 4 minutes
- User selectable languages
- Stores up to 7,000 results
- Barcode reader to scan calibration data, patient and operator ID
- Optional label printer available
- Optimised for Hospitals, Doctors and Pharmacies

Quo-Lab

▼ Product video:



¹⁾ Distributed by IME-DC GmbH

INSULIN PUMP



DANA Diabecare RS Smart Insulin Pump³⁾

- Dimension of Insulin Pump only 92 x 45 x 19 mm
- Weight of Insulin Pump only 62 g (including battery)
- Data Transfer via Bluetooth Low Energy (BLE)
- Access all pump functions from smartphone (Android and iOS applications)
- Extended bolus (30 minutes – 8 hours)
- Insulin reservoir capacity up to 300 units
- Water Protection Rating IP28 (protected from undefined immersion)
- Carbo / Bolus calculation programme
- Bolus from as small as 0.05 unit
- Minimum Basal 0.04 unit/hour
- Colours available in white, anthracite, black, green or pink

▼ Product video:



BLUETOOTH LOW ENERGY

AnyDANA - A
android

AnyDANA - i
iOS

Free app download to operate the insulin pump:

- **AnyDANA-A:** Google Play Store (minimum OS requirement: Android 4.3)
- **AnyDANA-i:** iTunes Store (minimum OS requirement: iOS 9.0)

Note: The pictured smartphones are not part of the sales product.

³⁾ Distributed by IME-DC GmbH

INSULIN THERAPY – INFUSION SETS

Easy-Release (Steel Cannula)⁴⁾

- Slide mechanism: Easily to connect and release
- Puncture angle: 90°
- Cannula: Stainless steel (27G)
- Self-adhesive patch: Enables the firm fixation with the puncture site
- *universal* – Compatible with most insulin pumps with a Luer port

The logo for SOOIL, featuring the word "SOOIL" in a blue, sans-serif font with a stylized blue starburst graphic above the letter "O".

Soft-Release-O (Teflon Cannula)⁴⁾

- Clip mechanism: Easily to connect and release
(directly on the cannula body / 360° rotatable)
- Puncture angle: 90°
- Cannula: Polyurethane (24G)
- Self-adhesive anti-allergy patch: Enables the firm fixation with the puncture site
- Suitable for all **DANA** Insulin Pumps

For further cannulas please visit our website.



Easy-Release (Steel Cannula)



Soft-Release-O (Teflon Cannula)

INSULIN THERAPY – PEN NEEDLES

IME-FINE

NEW

IME-FINE Universal Pen Needles

- Available sizes: 31G, 4 mm | 31G, 5 mm | 31G, 6 mm | 31G, 8 mm | 29G, 12 mm
- „Soft feel“-Cut for a soft insertion
- Silicone Coating for a gentle slide into the skin
- “Thin Wall”-Technology for a maximum insulin flow due to an increased inside diameter
- Available packaging unit: 100 pcs.

31G



Note: The pictured insulin pen is not part of the sales product.


INSULIN THERAPY – PEN NEEDLES

DEXFINE Universal Pen Needles

- Available sizes: 32G, 4 mm | 32G, 5 mm | 32G, 6 mm | 31G, 8 mm
- „Soft feel“-Cut for a soft insertion
- Silicone Coating for a gentle slide into the skin
- “Thin Wall”-Technology for a maximum insulin flow due to an increased inside diameter
- Available packaging unit: 100 pcs.



32G



LANCING DEVICE

IME-Sensitive

IME-Sensitive Lancing Device

- UNIVERSAL – Compatible with the most blood lancets*)
- Five variable depth of puncture for your personal preference
- Ergonomic design for a better handling
- One-hand usage possible
- Exact leading of lancets for smoother insertion
- Lancet-Ejection-Function for your safety



*) Compatible with blood lancets with outer diameter of the lancet body: 6.2 mm - 6.5 mm

IME-DC Blood Lancets

- UNIVERSAL – Compatible with the most lancing devices^{**})
- Gentle Penetration for a lower sensation of pain
- Easy Handling facilitates the process of blood collection
- Superfine Lancet 30G (0.30 mm)
- Protective Cap for a secure disposal
- Available packaging units: 50 pcs. | 100 pcs.

IME-DC

NEW



^{**}) Compatible with lancing devices with inner diameter of the lancet holder: 5.8 mm - 6.3 mm

SAFETY PRODUCT



DEXFINE protect Safety Pen Needles⁵⁾

- Visible Pen Needle Tip enables easy checking whether the pen is ready for injection
- Automatic Interlock immediately after injection
- User Protection against stab wounds
- Available sizes: 30G, 5 mm | 30G, 8 mm
- Available packaging unit: 100 pcs.



**COMING SOON
IN 2020**

⁵⁾ Trademark owned and distributed by IME-DC GmbH

DEXLANCET protect Safety Single Use Lancets⁶⁾

- Extraction of capillary blood by contact-activation
- Available size: 28G, 1.8 mm (other sizes on request)
- Automatic Interlock immediately after pricking
- User Protection against stab wounds
- Easy Use
- Available packaging unit: 100 pcs.



⁶⁾ Trademark owned and distributed by IME-DC GmbH

SKIN CARE PRODUCT



DiaAngel Cream Mousse

- Daily care for sensitive, dry and cracked skin
- Suitable for people with diabetes
- With 10 % Urea for an intensive moisturizing care
- Free of aromatic substances
- Reduces itchiness and prevents ulcerations
- Dermatologically tested
- Care effect has been clinically proven
- Filling quantity: 125 mL per can
- Made in Germany



iAngel Cream Mousse

- For tired and heavy legs
- Contains extracts of red grape leaf, horse chestnut and menthol
- With 3 % Urea for a moisturizing care
- Dermatologically tested
- Filling quantity: 125 mL per can
- Made in Germany

iAngel



DISINFECTION & HYGIENIC DISPOSAL



IME Swaps

- Cleaning Wipes for the skin
- Disinfection before needle injection
- With 70 % Isopropanol
- Available packaging unit: 100 pcs.

Disposal Container⁷⁾

- Up to 0.43 L volume
 - Round opening
 - Removal and rotation aid
 - In accordance with TRBA 250, UN, DIN EN ISO 23907 *
-
- Up to 1.53 L volume
 - Oval opening
 - Removal and rotation aid
 - In accordance with TRBA 250, AFNOR, UN, DIN EN ISO 23907 *



⁷⁾ Distributed by IME-DC GmbH | * Information supplied by the manufacturer (Sarstedt | Status: 2018-10)

SOFTWARE SOLUTION

DIABASS Management programme⁸⁾

- Diabetes diary for blood glucose, insulin or sporting activities
- Suitable for iDia and DEXTRA devices
- Transfer your blood glucose measurements to PC by USB data cable
- Visualisation of your blood glucose measurements by diagrams and charts

iDia

DEXTRA



⁸⁾ Product by mediaspects GmbH | <https://www.diabass.com/>

WHAT WE ARE LOOKING FOR

In view of our planned future growth and our goal of product availability in more than 50 countries, we are looking for professional distributors in the following regions:

- Europe
- Middle and South America
- Middle East and Northern Africa
- Africa
- Asia


You should meet the following requirements:



- Ability to implement country-specific business strategies



- Existing sales network with pharmacies and pharmacy wholesalers for guaranteed product availability



- Existing distribution network with doctors and hospitals specialised in diabetes therapy



- Existing sales team sufficient for achieving business strategy



- Financial stability necessary to realise projects



- High motivation to obtain significant market shares



- Experience with tender business

WHAT WE OFFER



- Portfolio of high quality products from a single source



- Patented new technologies and high quality solutions



- Exclusive distribution possible



- Co-Branding services for our customers



- Competitive prices




- Business strategies for market success



- High profit possibilities (country-specific)



- Partial investment support to realise common goals



- Service support directly from Germany



- Education and training courses for your team (product expertise as well as sales)



For further questions regarding a cooperation check, please contact:

Matthias Stark

Head of Sales International

E-mail: matthias.stark@imedc.de

Phone: +49 9281 | 85 01 6-0



REFERENCES

Partnerships



- Business cooperation in over 35 countries



- Cooperation with all German health insurance providers



- Contractual partner of health care organisations, hospitals, outpatient care and retirement homes



- Technology partner with Samsung SDS

Studies



- Interference studies



- Clinical study reports

Certifications




- EN ISO 13485

Qualifications



- PQS
Prequalification certificate for German Health Care System



- Work & Health Protection
Concept for employee's health, safety and work protection

Awards



- Mittelstandslöwe
Regional honour for successful companies in 2008



- CrefoZert
Solvency Certificate
www.creditreform.de/crefozert



- Innovationspreis-IT
E-health project „move...“
honoured with predicate „Best of 2016“ innovation




- Innovationspreis-IT
E-health project „move...“
honoured with predicate „Best of 2017“ innovation

For further questions regarding a cooperation check,
please do not hesitate to contact us!

For more information please visit our website
www.imedc.de

or scan the QR code with your mobile device.



IME-DC GmbH
Fuhrmannstraße 11
95030 Hof | Germany

Tel: +49 9281 | 85 01 6-0
Fax: +49 9281 | 85 01 6-100
E-Mail: info@imedc.de
Web: www.imedc.de